

Metafectene technical note for lentiviral production

Adriana Grbenicek¹, Bernhard Berkus², Patrick Maier¹

¹Department of Radiation Oncology, University Medical Center Mannheim, Heidelberg University, Germany

²DKFZ (German Cancer Research Centre), Heidelberg, Germany

Virus production: Transfection of 293T/17 cells for production of lentiviruses

The here described procedure for production of lentiviruses extends over four days plus three more days for determining the viral titer as described in [1].

Material:

Plasmid pHR'SINcPPT-SEW carrying EGFP as transgene controlled by the SFFV promoter Plasmid pCMVΔR8.91 packaging plasmid containing the HIV genes gag and pol Plasmid pMD.G encoding the envelope protein VSV-G

Metafectene	Biontex T020-5	
Na-butyrat	Sigma Aldrich B-5887	1 M stock solution in PBS
Poly-D-Lysine	Sigma P-7280	1 mg/ml stock solution in PBS
DMEM with FCS and P/S		
DMEM without additives		
Millex-Ha Filter 0.45 µm	Millipore SLHA 033SS	
Vivaspin 100 000 MW cut-off	Sartorius VS2041	

The lentiviral vector system

The lentiviral plasmid pHR'SINcPPT-SEW [2] contains the information of the recombinant lentiviral genome which is devoid of any information for the original HIV proteins. Instead, flanked by the 5'LTR (long terminal repeat domain) and the Δ 3'LTR (with deleted promoter sequences in order to reduce the risk of activation of downstream genes after integration into the human genome) the transgene EGFP (enhanced green fluorescent protein) is placed under the control of the U3 region of the LTR of SFFV (spleen focus-forming virus) and upstream of the WPRE (woodchuck hepatitis virus posttranscriptional regulatory element). An additional element is the packaging signal Ψ which is required that the viral RNA will be used as viral genome and packed into the viral capsid. The capsid protein and the reverse transcriptase are encoded by the genes gag and pol on the packaging plasmid pCMV Δ R8.91. The third plasmid delivers the information for the envelope protein VSV-G (vesicular stomatitis virus glycoprotein) for pseudotyping (enlargement of the spectrum of transducable cell types even of different species) of the viral capsid. This is a packaging system of the second generation (for further information about lentiviral vector systems see [3]).

Protocol

All information is for a culture dish of a diameter of 10 cm, and all solutions should be warmed up to room temperature.



Day 1: The culture dish is coated with 3 ml Poly-D-Lysine (0.1 mg/ml in PBS). After 5 min, Poly-D-Lysine is removed and stored at 4 °C since it can be used several times. The dish is washed with 5 ml sterile water and thus ready for use. 5 x 10^6 293T/17 cells are seeded in 14 ml DMEM with additives and cultured overnight.

Day 2: On the second day transfection is performed with three plasmids.

Overall 10 μ g of DNA is used for one culture dish: 4.4 μ g lentiviral plasmid, 3.4 μ g pCMV Δ R8.91, and 2.2 μ g pMD.G. The plasmids are given into a well of a 24-wellplate and mixed thoroughly with 700 μ l DMEM without additives \rightarrow solution A.

46 µl **Metafectene** is mixed in another well of a 24-wellplate with 700 µl of DMEM without additives \rightarrow solution B

Then solution A and B are mixed in one well and incubated for 20 min at room temperature. The solution of nearly 1.5 ml is added drop by drop to the cells with concomitant gentle orbital shaking of the dish. This is followed by an overnight incubation at 37 °C.

Day 3: In the morning, the medium is replaced by 14 ml DMEM with additives and 140 μ l Na-butyrate (10 mM).

In the evening the medium is replaced by 7 ml DMEM without additives.

Day 4: The medium is collected (it can be stored at 4 °C for several hours) and filtered with a 0.45 μ m filter. The cleared virus-containing solution is put on a Vivaspin-tube and concentrated up to 0.5 ml or 1 ml final volume (depending on the number of 10 cm dishes used for virus production) with 3000 rpm (\triangleq ca. 1800g) at 4 °C in an Eppendorf tabletop centrifuge. Aliquots of 50 μ l or 100 μ l are stored at -70 °C.



Figure 1. EGFP-expressing cells two days after transfection. Overlay (C) of the fluorescent channel (A) with the bright field picture (B) illustrates that a nearly 100 % transfection efficiency was achieved with the transfection conditions as described in our protocol.

Virus titration: Transduction of HT1080 cells

Material:24 well platepolybrene hexadimethrine bromideSigma H-9268polybrene working solution 800 µg/mlPBS sterile filtrated, storage at 4°C for monthsPBS without Ca/MgGibco 14190-094DMEMGibco 41965-039

Protocol All work must be done under S2 conditions.



Day 1: 8×10^4 HT 1080 cells are seeded per well of a24 well plate. Titration of each dilution should be done in duplicates or triplicates!

Day 2: Four dilutions are analysed: 1:100, 1:1000, 1:10000, 1:100000:

1. dilution (1:100): 10 μl viral supernatant + 980 μl medium + 10 μl Polybrene working solution

further dilutions: 100 μl of the previous dilution $\,+\,890$ μl Medium + 10 μl Polybrene working solution

The medium on cells is replaced by 300 μ l of each dilution.

Day 3: The medium is exchanged with 0.5 ml DMEM without Polybrene per well

Day 4 or 5: The cells are harvested, washed once with 1 ml PBS (centrifuge 3 min with

3000 rpm \triangleq ca. 1800g in an Eppendorf tabletop centrifuge) and 0.3-1 x10⁶ cells are resuspended in 300 µl staining medium for FACS.

Analysis

Example: 80000 cells $\triangleq 100 \%$

FACS result: 22.5 % GFP-positive cells 22.5 % transduction efficiency

supernatant diluted 1:1000 in 300 μl

=> amount of viral particles in 1 ml concentrated supernatant = (transduction efficiency x total cell number)/100 % x (1 ml/volume of medium containing diluted viral particles) x dilution factor

= $(22.5 \% x 80000) / 100 \% x (1000 \mu l/300 \mu l) x 1000$

= 22.5 x 8000/3 x 1000

 $= 6 \text{ x} 10^7$ infectious particles (titer) per 1 ml

References:

- 1. Maier, P.; Herskind, C.; Fleckenstein, K.; Spier, I.; Laufs, S.; Zeller, W. J.; Fruehauf, S.; Wenz, F., MDR1 gene transfer using a lentiviral SIN vector confers radioprotection to human CD34+ hematopoietic progenitor cells. *Radiat Res* 2008, 169, (3), 301-10.
- 2. Demaison, C.; Parsley, K.; Brouns, G.; Scherr, M.; Battmer, K.; Kinnon, C.; Grez, M.; Thrasher, A. J., High-level transduction and gene expression in hematopoietic repopulating cells using a human immunodeficiency [correction of imunodeficiency] virus type 1-based lentiviral vector containing an internal spleen focus forming virus promoter. *Hum Gene Ther* **2002**, 13, (7), 803-13.
- 3. Maier, P.; von Kalle, C.; Laufs, S., Retroviral vectors for gene therapy. *Future Microbiol* **2010**, 5, 1507-23.